



IMPORTANT NOTICE

WICHTIGER HINWEIS / REMARQUE IMPORTANTE
 NOTA IMPORTANTE / NOTA IMPORTANTE
 NOTA IMPORTANTE / VIKTIG ANMÄRKNING
 VIGTIGT / DŮLEŽITÉ UPOZORNĚNÍ
 DŮLEŽITÉ UPOZORNENIE / WAŻNE INFORMACJE
 FONTOS MEGJEGYZÉS / ΣΗΜΑΝΤΙΚΗ ΣΗΜΕΙΩΣΗ

Bilirubin Total DPD Gen.2 COBAS INTEGRA 400/700/800

03146022 322, Sys.-ID 07 6587 2

BILT2: Test-ID 0-587

Please refer to the new assigned values:
 Bitte beachten Sie die neu ermittelten Sollwerte:
 Veuillez tenir compte des nouvelles valeurs théoriques:
 Observe los nuevos valores teóricos:
 Si prega di prendere in considerazione i nuovi valori teorici:
 Consulte os novos valores teóricos:
 Vänligen beakta de nya åsatta värdena:
 Se de nye angivne værdier:
 Věnujte pozornost novým určeným hodnotám :
 Prosím použite novo definované hodnoty:
 Należy odnieść się do nowo wyznaczonych wartości:
 Nézze meg az új vonatkozó értékeket:
 Ανατρέξτε στις νέες καθορισμένες τιμές:

Plasma / Plazma / Osocze / Πλάσμα
 Serum / SÉRUM / Suero / Siero / Soro / SÉRUM / Surowica / Szérum / Ορός

LOT	μmol/L	mg/dL	mg/L
	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s., No. 10759350 190

168934	78.8	4.61	46.1
168935	82.2	4.81	48.1
168936	83.8	4.90	49.0
169826	83.2	4.87	48.7
169827	81.0	4.74	47.4
169828	82.0	4.80	48.0
169829	81.2	4.75	47.5
169830	78.5	4.59	45.9
169831	80.2	4.69	46.9
170061	87.4	5.11	51.1
170062	86.1	5.04	50.4
170063	86.8	5.08	50.8
170064	84.4	4.94	49.4
170065	85.7	5.01	50.1
170066	83.7	4.90	49.0
171683	85.8	5.02	50.2
171684	85.0	4.97	49.7
171685	82.7	4.84	48.4
171686	82.7	4.84	48.4
171687	83.3	4.87	48.7
171688	84.3	4.93	49.3

LOT	µmol/L	mg/dL	mg/L
	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvårde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvårde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvårde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s., No. 10759350 190

172725	88.3	5.17	51.7
172726	85.2	4.98	49.8
172727	87.0	5.09	50.9
172728	87.0	5.09	50.9
172729	86.5	5.06	50.6
172730	86.2	5.04	50.4
173630	78.1	4.57	45.7
173631	76.6	4.48	44.8
173632	79.2	4.63	46.3
173633	78.6	4.60	46.0
173634	78.6	4.60	46.0
173635	77.7	4.55	45.5
173928	96.5	5.65	56.5
173929	95.4	5.58	55.8

LOT	Roche µmol/L			Roche mg/dL			Roche mg/L		
	Value/Wert Valeur/Valor Valore/Valor Vårde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD	Value/Wert Valeur/Valor Valore/Valor Vårde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD	Value/Wert Valeur/Valor Valore/Valor Vårde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD

Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)

167326	16.2	13.2 - 19.2	1.0	0.948	0.777 - 1.119	0.057	9.48	7.77 - 11.19	0.57
167486	17.6	14.3 - 20.9	1.1	1.03	0.85 - 1.21	0.06	10.3	8.5 - 12.1	0.6
167969	18.3	15.0 - 21.6	1.1	1.07	0.89 - 1.25	0.06	10.7	8.9 - 12.5	0.6
167972	17.8	14.5 - 21.1	1.1	1.04	0.86 - 1.22	0.06	10.4	8.6 - 12.2	0.6
168314	17.6	14.3 - 20.9	1.1	1.03	0.85 - 1.21	0.06	10.3	8.5 - 12.1	0.6
168317	17.5	14.2 - 20.8	1.1	1.02	0.84 - 1.20	0.06	10.2	8.4 - 12.0	0.6
170571	19.4	15.8 - 23.0	1.2	1.13	0.92 - 1.34	0.07	11.3	9.2 - 13.4	0.7
170574	18.6	15.3 - 21.9	1.1	1.09	0.88 - 1.30	0.07	10.9	8.8 - 13.0	0.7
171074	18.6	15.3 - 21.9	1.1	1.09	0.88 - 1.30	0.07	10.9	8.8 - 13.0	0.7
171075	19.4	15.8 - 23.0	1.2	1.13	0.92 - 1.34	0.07	11.3	9.2 - 13.4	0.7
172101	18.8	15.5 - 22.1	1.1	1.10	0.89 - 1.31	0.07	11.0	8.9 - 13.1	0.7
172104	19.5	15.9 - 23.1	1.2	1.14	0.93 - 1.35	0.07	11.4	9.3 - 13.5	0.7
173499	18.2	14.9 - 21.5	1.1	1.06	0.88 - 1.24	0.06	10.6	8.8 - 12.4	0.6
173501	17.8	14.5 - 21.1	1.1	1.04	0.86 - 1.22	0.06	10.4	8.6 - 12.2	0.6
175650	18.5	15.2 - 21.8	1.1	1.08	0.90 - 1.26	0.06	10.8	9.0 - 12.6	0.6
176136	18.0	14.7 - 21.3	1.1	1.05	0.87 - 1.23	0.06	10.5	8.7 - 12.3	0.6
176138	18.0	14.7 - 21.3	1.1	1.05	0.87 - 1.23	0.06	10.5	8.7 - 12.3	0.6
176469	17.6	14.3 - 20.9	1.1	1.03	0.85 - 1.21	0.06	10.3	8.5 - 12.1	0.6

LOT	Roche µmol/L			Roche mg/dL			Roche mg/L		
	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD
	Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE
	Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP
	Värde/Værdi	Område/Område	1DP/1SD	Värde/Værdi	Område/Område	1DP/1SD	Värde/Værdi	Område/Område	1DP/1SD
	Homota/Hodnota	Rozmezí/Rozsah	1SD/1SD	Homota/Hodnota	Rozmezí/Rozsah	1SD/1SD	Homota/Hodnota	Rozmezí/Rozsah	1SD/1SD
	Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD
	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

Precinorm U plus, No. 12149435 122

168539	15.7	13.0 - 18.4	0.9	0.918	0.753 - 1.083	0.055	9.18	7.53 - 10.83	0.55
172222	16.3	13.3 - 19.3	1.0	0.954	0.783 - 1.125	0.057	9.54	7.83 - 11.25	0.57
175487	17.9	14.6 - 21.2	1.1	1.05	0.87 - 1.23	0.06	10.5	8.7 - 12.3	0.6

Precipath U, No. 10171760 122, 10171778 122, 10651265 922 (Q.C.S.)

167436	62.8	51.4 - 74.2	3.8	3.67	3.01 - 4.33	0.22	36.7	30.1 - 43.3	2.2
168009	65.3	53.6 - 77.0	3.9	3.82	3.13 - 4.51	0.23	38.2	31.3 - 45.1	2.3
168309	69.9	57.3 - 82.5	4.2	4.09	3.34 - 4.84	0.25	40.9	33.4 - 48.4	2.5
169839	85.2	69.9 - 100.5	5.1	4.98	4.08 - 5.88	0.30	49.8	40.8 - 58.8	3.0
169840	87.2	71.6 - 102.8	5.2	5.10	4.17 - 6.03	0.31	51.0	41.7 - 60.3	3.1
171368	84.4	69.1 - 99.7	5.1	4.94	4.04 - 5.84	0.30	49.4	40.4 - 58.4	3.0
173043	91.0	74.5 - 107.5	5.5	5.32	4.36 - 6.28	0.32	53.2	43.6 - 62.8	3.2
174531	90.3	74.1 - 106.5	5.4	5.28	4.32 - 6.24	0.32	52.8	43.2 - 62.4	3.2
175332	81.7	67.0 - 96.4	4.9	4.78	3.91 - 5.65	0.29	47.8	39.1 - 56.5	2.9
175333	83.1	68.1 - 98.1	5.0	4.86	3.99 - 5.73	0.29	48.6	39.9 - 57.3	2.9
176287	78.0	63.9 - 92.1	4.7	4.56	3.75 - 5.37	0.27	45.6	37.5 - 53.7	2.7

Precipath U plus, No. 12149443 122

166448	79.4	65.0 - 93.8	4.8	4.64	3.80 - 5.48	0.28	46.4	38.0 - 54.8	2.8
171737	79.0	64.9 - 93.1	4.7	4.62	3.78 - 5.46	0.28	46.2	37.8 - 54.6	2.8
176975	81.3	66.6 - 96.0	4.9	4.76	3.89 - 5.63	0.29	47.6	38.9 - 56.3	2.9
176976	82.6	67.6 - 97.6	5.0	4.83	3.96 - 5.70	0.29	48.3	39.6 - 57.0	2.9
176977	82.3	67.6 - 97.0	4.9	4.81	3.94 - 5.68	0.29	48.1	39.4 - 56.8	2.9

Precibil, No. 10158046 122

170941	218	179 - 257	13	12.8	10.4 - 15.2	0.8	128	104 - 152	8
177006	226	184 - 268	14	13.2	10.8 - 15.6	0.8	132	108 - 156	8

LOT	µmol/L		mg/dL		mg/L	
	Sollwert	Bereich BÄK	Sollwert	Bereich BÄK	Sollwert	Bereich BÄK

Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)

167326	16.2	9.4 - 23.0	0.948	0.548 - 1.348	9.48	5.48 - 13.48
167486	17.6	10.8 - 24.4	1.03	0.63 - 1.43	10.3	6.3 - 14.3
167969	18.3	11.5 - 25.1	1.07	0.67 - 1.47	10.7	6.7 - 14.7
167972	17.8	11.0 - 24.6	1.04	0.64 - 1.44	10.4	6.4 - 14.4
168314	17.6	10.8 - 24.4	1.03	0.63 - 1.43	10.3	6.3 - 14.3
168317	17.5	10.7 - 24.3	1.02	0.62 - 1.42	10.2	6.2 - 14.2
170571	19.4	12.6 - 26.2	1.13	0.73 - 1.53	11.3	7.3 - 15.3
170574	18.6	11.8 - 25.4	1.09	0.69 - 1.49	10.9	6.9 - 14.9
171074	18.6	11.8 - 25.4	1.09	0.69 - 1.49	10.9	6.9 - 14.9
171075	19.4	12.6 - 26.2	1.13	0.73 - 1.53	11.3	7.3 - 15.3
172101	18.8	12.0 - 25.6	1.10	0.70 - 1.50	11.0	7.0 - 15.0
172104	19.5	12.7 - 26.3	1.14	0.74 - 1.54	11.4	7.4 - 15.4

LOT	µmol/L				mg/dL				mg/L			
	Sollwert	Bereich BÄK			Sollwert	Bereich BÄK			Sollwert	Bereich BÄK		

Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)

173499	18.2	11.4	-	25.0	1.06	0.66	-	1.46	10.6	6.6	-	14.6
173501	17.8	11.0	-	24.6	1.04	0.64	-	1.44	10.4	6.4	-	14.4
175650	18.5	11.7	-	25.3	1.08	0.68	-	1.48	10.8	6.8	-	14.8
176136	18.0	11.2	-	24.8	1.05	0.65	-	1.45	10.5	6.5	-	14.5
176138	18.0	11.2	-	24.8	1.05	0.65	-	1.45	10.5	6.5	-	14.5
176469	17.6	10.8	-	24.4	1.03	0.63	-	1.43	10.3	6.3	-	14.3

Precinorm U plus, No. 12149435 122

168539	15.7	8.9	-	22.5	0.918	0.518	-	1.318	9.18	5.18	-	13.18
172222	16.3	9.5	-	23.1	0.954	0.554	-	1.354	9.54	5.54	-	13.54
175487	17.9	11.1	-	24.7	1.05	0.65	-	1.45	10.5	6.5	-	14.5

Precipath U, No. 10171760 122, 10171778 122, 10651265 922 (Q.C.S.)

167436	62.8	46.6	-	79.0	3.67	2.71	-	4.63	36.7	27.1	-	46.3
168009	65.3	48.2	-	82.4	3.82	2.83	-	4.81	38.2	28.3	-	48.1
168309	69.9	51.6	-	88.2	4.09	3.04	-	5.14	40.9	30.4	-	51.4
169839	85.2	63.0	-	107.4	4.98	3.69	-	6.27	49.8	36.9	-	62.7
169840	87.2	64.4	-	110.0	5.10	3.78	-	6.42	51.0	37.8	-	64.2
171368	84.4	62.5	-	106.3	4.94	3.65	-	6.23	49.4	36.5	-	62.3
173043	91.0	67.3	-	114.7	5.32	3.94	-	6.70	53.2	39.4	-	67.0
174531	90.3	66.9	-	113.7	5.28	3.90	-	6.66	52.8	39.0	-	66.6
175332	81.7	60.4	-	103.0	4.78	3.55	-	6.01	47.8	35.5	-	60.1
175333	83.1	61.5	-	104.7	4.86	3.60	-	6.12	48.6	36.0	-	61.2
176287	78.0	57.6	-	98.4	4.56	3.36	-	5.76	45.6	33.6	-	57.6

Precipath U plus, No. 12149443 122

166448	79.4	58.7	-	100.1	4.64	3.44	-	5.84	46.4	34.4	-	58.4
171737	79.0	58.6	-	99.4	4.62	3.42	-	5.82	46.2	34.2	-	58.2
176975	81.3	60.3	-	102.3	4.76	3.53	-	5.99	47.6	35.3	-	59.9
176976	82.6	61.0	-	104.2	4.83	3.57	-	6.09	48.3	35.7	-	60.9
176977	82.3	61.0	-	103.6	4.81	3.55	-	6.07	48.1	35.5	-	60.7

Precibil, No. 10158046 122

170941	218	161	-	275	12.8	9.5	-	16.1	128	95	-	161
177006	226	166	-	286	13.2	9.9	-	16.5	132	99	-	165