

CHE2/CHED2

Note 2006-03- 04498577 / 04498631 - 1 - 1



IMPORTANT NOTICE

WICHTIGER HINWEIS / REMARQUE IMPORTANTE
NOTA IMPORTANTE / NOTA IMPORTANTE
NOTA IMPORTANTE / VIKTIG ANMÄRKNING
VIGTIGT / DŮLEŽITÉ UPOZORNĚNÍ
DŮLEŽITÉ UPOZORNENIE / WAŽNE INFORMACJE
FONTOS MEGJEGYZÉS / ΣΗΜΑΝΤΙΚΗ ΣΗΜΕΙΩΣΗ

Please refer to the new assigned values:
Bitte beachten Sie die neu ermittelten Sollwerte:
Veuillez tenir compte des nouvelles valeurs théoriques:
Observe los nuevos valores teóricos:
Si prega di prendere in considerazione i nuovi valori teorici:
Consulte os novos valores teóricos:
Vänligen beakta de nya åsatta värdena:
Se de nye angivne værdier:
Věnujte pozornost novým určeným hodnotám :
Prosím použite novo definované hodnoty:
Należy odnieść się do nowo wyznaczonych wartości:
Nézze meg az új vonatkozó értékeket:
Ανατρέξτε στις νέες καθορισμένες τιμές:

Cholinesterase Gen.2

04498577 190, Sys.-ID 07 6842 1

Cholinesterase/Dibucaine Gen.2

04498631 190, Sys.-ID 07 6845 6

COBAS INTEGRA 400/700/800

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός
Plasma / Plazma / Osocze / Πλάσμα

LOT	U/L	µkat/L
	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärderi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärderi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s., No. 10759350 190

167216	4180	69.8
167217	4230	70.6
167218	4200	70.1
167219	4200	70.1
167220	4210	70.3
167221	4200	70.1
168934	4430	74.0
168935	4420	73.8
168936	4440	74.1
169826	3850	64.3
169827	3860	64.5
169828	3860	64.5
169829	3840	64.1
169830	3860	64.5
169831	3830	64.0

CHE2/CHED2

Note 2006-03- 04498577 / 04498631 - 2 - 1

LOT	U/L	µkat/L
	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsværdi/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsværdi/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s., No. 10759350 190

170061	5110	85.3
170062	5080	84.8
170063	5090	85.0
170064	5080	84.8
170065	5080	84.8
170066	5050	84.3
171683	4820	80.5
171684	4830	80.7
171685	4800	80.2
171686	4870	81.3
171687	4840	80.8
171688	4850	81.0

LOT	Roche U/L			Roche µkat/L		
	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD
Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE	
Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP	
Værde/Værdi	Område/Område	1DP/1SD	Værde/Værdi	Område/Område	1DP/1SD	
Hodnota/Hodnota	Rozmezí/Rozsah	1SD/1SD	Hodnota/Hodnota	Rozmezí/Rozsah	1SD/1SD	
Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD	
Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	

Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)

167326	5320	4360 - 6280	320	88.8	72.9 - 104.7	5.3
167486	4480	3670 - 5290	270	74.8	61.3 - 88.3	4.5
167969	4880	4010 - 5750	290	81.5	66.8 - 96.2	4.9
167972	4880	4010 - 5750	290	81.5	66.8 - 96.2	4.9
168314	5730	4710 - 6750	340	95.7	78.6 - 112.8	5.7
168317	5730	4710 - 6750	340	95.7	78.6 - 112.8	5.7
170571	5980	4900 - 7060	360	99.9	81.9 - 117.9	6.0
170574	5980	4900 - 7060	360	99.9	81.9 - 117.9	6.0
171074	5910	4860 - 6960	350	98.7	81.0 - 116.4	5.9
171075	5910	4860 - 6960	350	98.7	81.0 - 116.4	5.9
171680	6260	5120 - 7400	380	105	87 - 123	6
172101	5440	4450 - 6430	330	90.8	74.6 - 107.0	5.4
172104	5440	4450 - 6430	330	90.8	74.6 - 107.0	5.4
173499	4180	3430 - 4930	250	69.8	57.2 - 82.4	4.2
173501	4180	3430 - 4930	250	69.8	57.2 - 82.4	4.2

CHE2/CHED2

Note 2006-03-04498577 / 04498631 - 3 - 1

LOT	Roche U/L			Roche µkat/L		
	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD
	Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE
	Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP
	Värde/Værdi	Område/Område	1DP/1SD	Värde/Værdi	Område/Område	1DP/1SD
	Hornota/Hodnota	Rozmezí/Rozsah	1SD/1SD	Hornota/Hodnota	Rozmezí/Rozsah	1SD/1SD
	Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD
	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

Precinorm U plus, 12149435 122

166002	5520	4530 - 6510	330	92.2	75.7 - 108.7	5.5
168539	4110	3360 - 4860	250	68.6	56.3 - 80.9	4.1
172222	4790	3920 - 5660	290	80.0	65.6 - 94.4	4.8

Precipath U, No. 10171760 122, 10171778 122, 10651265 922 (Q.C.S.)

167436	4650	3810 - 5490	280	77.7	63.6 - 91.8	4.7
168009	5000	4100 - 5900	300	83.5	68.5 - 98.5	5.0
168309	5400	4440 - 6360	320	90.2	74.0 - 106.4	5.4
169210	5310	4350 - 6270	320	88.7	72.8 - 104.6	5.3
169839	5690	4670 - 6710	340	95.0	77.9 - 112.1	5.7
169840	5690	4670 - 6710	340	95.0	77.9 - 112.1	5.7
171368	4820	3950 - 5690	290	80.5	66.1 - 94.9	4.8
171681	5710	4690 - 6730	340	95.4	78.3 - 112.5	5.7
173043	3640	2980 - 4300	220	60.8	50.0 - 71.6	3.6

Precipath U plus, 12149443 122

166448	5050	4150 - 5950	300	84.3	69.0 - 99.6	5.1
171737	4600	3760 - 5440	280	76.8	63.0 - 90.6	4.6

LOT	U/L		µkat/L	
	Sollwert	Bereich BÄK	Sollwert	Bereich BÄK

Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)

167326	5320	4360 - 6280	88.8	72.9 - 104.7
167486	4480	3670 - 5290	74.8	61.3 - 88.3
167969	4880	4010 - 5750	81.5	66.8 - 96.2
167972	4880	4010 - 5750	81.5	66.8 - 96.2
168314	5730	4710 - 6750	95.7	78.6 - 112.8
168317	5730	4710 - 6750	95.7	78.6 - 112.8
170571	5980	4900 - 7060	99.9	81.9 - 117.9
170574	5980	4900 - 7060	99.9	81.9 - 117.9
171074	5910	4860 - 6960	98.7	81.0 - 116.4
171075	5910	4860 - 6960	98.7	81.0 - 116.4
171680	6260	5120 - 7400	105	87 - 123
172101	5440	4450 - 6430	90.8	74.6 - 107.0
172104	5440	4450 - 6430	90.8	74.6 - 107.0
173499	4180	3430 - 4930	69.8	57.2 - 82.4
173501	4180	3430 - 4930	69.8	57.2 - 82.4

CHE2/CHED2

Note 2006-03- 04498577 / 04498631 - 4 - 1

LOT	U/L			µkat/L		
	Sollwert	Bereich BÄK		Sollwert	Bereich BÄK	

Precipath U, No. 10171760 122, 10171778 122, 10651265 922 (Q.C.S.)

167436	4650	3810	-	5490	77.7	63.6	-	91.8
168009	5000	4100	-	5900	83.5	68.5	-	98.5
168309	5400	4440	-	6360	90.2	74.0	-	106.4
169210	5310	4350	-	6270	88.7	72.8	-	104.6
169839	5690	4670	-	6710	95.0	77.9	-	112.1
169840	5690	4670	-	6710	95.0	77.9	-	112.1
171368	4820	3950	-	5690	80.5	66.1	-	94.9
171681	5710	4690	-	6730	95.4	78.3	-	112.5
173043	3640	2980	-	4300	60.8	50.0	-	71.6

COBAS INTEGRA, PRECINORM and PRECIPATH are trademarks of Roche.
©2006 Roche Diagnostics

2006-03



Roche Diagnostics GmbH, D-68298 Mannheim

For USA:
Roche Diagnostics, Indianapolis, IN
US Customer Support 1-800-428-2336

