


**IMPORTANT  
NOTICE**

**WICHTIGER HINWEIS / REMARQUE IMPORTANTE  
NOTA IMPORTANTE / NOTA IMPORTANTE  
NOTA IMPORTANTE / VIKTIG ANMÄRKNING  
VIGTIGT / DŮLEŽITÉ UPOZORNĚNÍ  
DŮLEŽITÉ UPOZORNENIE / WAŽNE INFORMACJE  
FONTOS MEGJEGYZÉS / ΣΗΜΑΝΤΙΚΗ ΣΗΜΕΙΩΣΗ**

**γ-Glutamyltransferase ver.2**  
Standardized against IFCC

**COBAS INTEGRA 400/400plus/700/800**

**03002721** 122, System-ID 07 6598 8

Please refer to the new assigned values:  
Bitte beachten Sie die neu ermittelten Sollwerte:  
Veuillez tenir compte des nouvelles valeurs théoriques:  
Observe los nuevos valores teóricos:  
Si prega di prendere in considerazione i nuovi valori teorici:  
Consulte os novos valores teóricos:  
Vänligen beakta de nya åsatta värdena:  
Se de nye angivne værdier:  
Věnujte pozornost novým určeným hodnotám :  
Prosím použite novo definované hodnoty:  
Należy odnieść się do nowo wyznaczonych wartości:  
Nézzze meg az új vonatkozó értékeket:  
Ανατρέξτε στις νέες καθορισμένες τιμές:

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός  
Plasma / Plazma / Osocze / Πλάσμα

LOT	U/L	µkat/L
	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvårde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvårde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s., No. 10759350 190

168934	126	2.10
168935	126	2.10
168936	126	2.10
169826	103	1.72
169827	102	1.70
169828	103	1.72
169829	103	1.72
169830	103	1.72
169831	103	1.72
170061	120	2.00
170062	121	2.02
170063	121	2.02
170064	121	2.02
170065	121	2.02
170066	120	2.00
171683	110	1.84
171684	110	1.84
171685	109	1.82
171686	109	1.82
171687	109	1.82
171688	110	1.84

LOT	Roche U/L			Roche µkat/L		
	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD
	Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE
	Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP
	Värde/Værdi	Område/Område	1DP/1SD	Värde/Værdi	Område/Område	1DP/1SD
	Hornota/Hodnota	Rozmezi/Rozsah	1SD/1SD	Hornota/Hodnota	Rozmezi/Rozsah	1SD/1SD
	Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD
	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

**Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)**

167326	46.7	38.3 - 55.1	2.8	0.780	0.639 - 0.921	0.047
167486	45.1	37.0 - 53.2	2.7	0.753	0.618 - 0.888	0.045
167969	46.3	37.9 - 54.7	2.8	0.773	0.635 - 0.911	0.046
167972	46.3	37.9 - 54.7	2.8	0.773	0.635 - 0.911	0.046
168314	46.2	37.8 - 54.6	2.8	0.772	0.634 - 0.910	0.046
168317	46.2	37.8 - 54.6	2.8	0.772	0.634 - 0.910	0.046
170571	50.8	41.8 - 59.8	3.0	0.848	0.695 - 1.001	0.051
170574	50.8	41.8 - 59.8	3.0	0.848	0.695 - 1.001	0.051
171074	49.3	40.3 - 58.3	3.0	0.823	0.676 - 0.970	0.049
171075	49.3	40.3 - 58.3	3.0	0.823	0.676 - 0.970	0.049
171680	45.5	37.4 - 53.6	2.7	0.760	0.622 - 0.898	0.046
172101	48.9	40.2 - 57.6	2.9	0.817	0.670 - 0.964	0.049
172104	48.9	40.2 - 57.6	2.9	0.817	0.670 - 0.964	0.049

**Precinorm U plus, 12149435 122**

166002	43.4	35.6 - 51.2	2.6	0.725	0.593 - 0.857	0.044
168539	44.3	36.2 - 52.4	2.7	0.740	0.608 - 0.872	0.044
172222	46.1	37.7 - 54.5	2.8	0.770	0.632 - 0.908	0.046

**Precipath U, No. 10171760 122, 10171778 122, 10651265 922 (Q.C.S.)**

167436	221	182 - 260	13	3.69	3.03 - 4.35	0.22
168009	223	184 - 262	13	3.72	3.06 - 4.38	0.22
168309	213	174 - 252	13	3.56	2.93 - 4.19	0.21
169210	209	170 - 248	13	3.49	2.86 - 4.12	0.21
169839	220	181 - 259	13	3.67	3.01 - 4.33	0.22
169840	220	181 - 259	13	3.67	3.01 - 4.33	0.22
171368	223	184 - 262	13	3.72	3.06 - 4.38	0.22
171681	208	172 - 244	12	3.47	2.84 - 4.10	0.21
173043	236	194 - 278	14	3.94	3.22 - 4.66	0.24

**Precipath U plus, 12149443 122**

166448	190	157 - 223	11	3.17	2.60 - 3.74	0.19
171737	210	171 - 249	13	3.51	2.88 - 4.14	0.21

# GGT12

Note 2006 - 04 - 03002721 - 3 - 3

LOT	U/L			µkat/L		
	Sollwert	Bereich BÄK		Sollwert	Bereich BÄK	

### Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)

167326	46.7	36.5	-	56.9	0.780	0.609	-	0.951
167486	45.1	35.2	-	55.0	0.753	0.588	-	0.918
167969	46.3	36.1	-	56.5	0.773	0.602	-	0.944
167972	46.3	36.1	-	56.5	0.773	0.602	-	0.944
168314	46.2	36.0	-	56.4	0.772	0.601	-	0.943
168317	46.2	36.0	-	56.4	0.772	0.601	-	0.943
170571	50.8	39.7	-	61.9	0.848	0.662	-	1.034
170574	50.8	39.7	-	61.9	0.848	0.662	-	1.034
171074	49.3	38.5	-	60.1	0.823	0.643	-	1.003
171075	49.3	38.5	-	60.1	0.823	0.643	-	1.003
171680	45.5	35.6	-	55.4	0.760	0.592	-	0.928
172101	48.9	38.1	-	59.7	0.817	0.637	-	0.997
172104	48.9	38.1	-	59.7	0.817	0.637	-	0.997

### Precipath U, No. 10171760 122, 10171778 122, 10651265 922 (Q.C.S.)

167436	221	173	-	269	3.69	2.88	-	4.50
168009	223	175	-	271	3.72	2.91	-	4.53
168309	213	165	-	261	3.56	2.78	-	4.34
169210	209	164	-	254	3.49	2.71	-	4.27
169839	220	172	-	268	3.67	2.86	-	4.48
169840	220	172	-	268	3.67	2.86	-	4.48
171368	223	175	-	271	3.72	2.91	-	4.53
171681	208	163	-	253	3.47	2.72	-	4.22
173043	236	185	-	287	3.94	3.07	-	4.81

PRECINORM, PRECIPATH and COBAS INTEGRA are trademarks of Roche.  
©2006 Roche Diagnostics

2006-04



Roche Diagnostics GmbH, D-68298 Mannheim

**For USA:**  
Roche Diagnostics, Indianapolis, IN  
US Customer Support 1-800-428-2336

