



## IMPORTANT NOTICE

**WICHTIGER HINWEIS / REMARQUE IMPORTANTE  
 NOTA IMPORTANTE / NOTA IMPORTANTE  
 NOTA IMPORTANTE / VIKTIG ANMÄRKNING  
 VIGTIGT / DŮLEŽITÉ UPOZORNĚNÍ  
 DŮLEŽITÉ UPOZORNENIE / WAŽNE INFORMACJE  
 FONTOS MEGJEGYZÉS / ΣΗΜΑΝΤΙΚΗ ΣΗΜΕΙΩΣΗ**

## Iron Gen.2 COBAS INTEGRA 400/700/800

03183696 122, Sys.-ID 07 6596 1

**IRON2** : Test-ID 0-596

Please refer to the new assigned values:  
 Bitte beachten Sie die neu ermittelten Sollwerte:  
 Veuillez tenir compte des nouvelles valeurs théoriques:  
 Observe los nuevos valores teóricos:  
 Si prega di prendere in considerazione i nuovi valori teorici:  
 Consulte os novos valores teóricos:  
 Vänligen beakta de nya åsatta värdena:  
 Se de nye angivne værdier:  
 Věnujte pozornost novým určeným hodnotám :  
 Prosím použite novo definované hodnoty:  
 Należy odnieść się do nowo wyznaczonych wartości:  
 Nézze meg az új vonatkozó értékeket:  
 Ανατρέξτε στις νέες καθορισμένες τιμές:

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός  
 Plasma / Plazma / Osocze / Πλάσμα

	<b>µmol/L</b>	<b>µg/dL</b>	<b>mg/L</b>
<b>LOT</b>	Calibration Value/Kalibrationswert	Calibration Value/Kalibrationswert	Calibration Value/Kalibrationswert
	Valeur de calibration/Valor de calibración	Valeur de calibration/Valor de calibración	Valeur de calibration/Valor de calibración
	Valore di calibrazione/Valor de calibração	Valore di calibrazione/Valor de calibração	Valore di calibrazione/Valor de calibração
	Kalibreringsvärde/Kalibreringsværdi	Kalibreringsvärde/Kalibreringsværdi	Kalibreringsvärde/Kalibreringsværdi
	Kalibrační hodnota/Hodnota kalibrácie	Kalibrační hodnota/Hodnota kalibrácie	Kalibrační hodnota/Hodnota kalibrácie
	Wartość kalibratora/Kalibrációs érték	Wartość kalibratora/Kalibrációs érték	Wartość kalibratora/Kalibrációs érték
	Τιμή βαθμονόμησης	Τιμή βαθμονόμησης	Τιμή βαθμονόμησης

**C.f.a.s., No. 10759350 190**

168934	33.6	188	1.88
168935	33.5	187	1.87
168936	33.4	187	1.87
169826	35.5	198	1.98
169827	35.5	198	1.98
169828	35.2	197	1.97
169829	35.6	199	1.99
169830	35.3	197	1.97
169831	35.3	197	1.97
170061	34.9	195	1.95
170062	34.7	194	1.94
170063	34.9	195	1.95
170064	34.9	195	1.95
170065	34.9	195	1.95
170066	34.8	195	1.95
171683	35.2	197	1.97
171684	35.2	197	1.97
171685	35.2	197	1.97
171686	35.2	197	1.97
171687	35.1	196	1.96
171688	35.3	197	1.97

LOT	<b>µmol/L</b>	<b>µg/dL</b>	<b>mg/L</b>
	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

**C.f.a.s., No. 10759350 190**

172725	35.8	200	2.00
172726	35.7	200	2.00
172727	35.8	200	2.00
172728	35.6	199	1.99
172729	35.5	198	1.98
172730	35.8	200	2.00
173630	36.0	201	2.01
173631	35.8	200	2.00
173632	35.8	200	2.00
173633	36.0	201	2.01
173634	35.8	200	2.00
173635	35.7	200	2.00
173928	33.5	187	1.87
173929	33.4	187	1.87

LOT	<b>Roche µmol/L</b>			<b>Roche µg/dL</b>			<b>Roche mg/L</b>		
	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD

**Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)**

167326	20.7	17.1 - 24.3	1.2	116	95 - 137	7	1.16	0.95 - 1.37	0.07
167486	20.5	16.9 - 24.1	1.2	115	94 - 136	7	1.15	0.94 - 1.36	0.07
167969	21.9	18.0 - 25.8	1.3	122	101 - 143	7	1.22	1.01 - 1.43	0.07
167972	21.9	18.0 - 25.8	1.3	122	101 - 143	7	1.22	1.01 - 1.43	0.07
168314	20.6	17.0 - 24.2	1.2	115	94 - 136	7	1.15	0.94 - 1.36	0.07
168317	20.6	17.0 - 24.2	1.2	115	94 - 136	7	1.15	0.94 - 1.36	0.07
170571	20.0	16.4 - 23.6	1.2	112	91 - 133	7	1.12	0.91 - 1.33	0.07
170574	20.0	16.4 - 23.6	1.2	112	91 - 133	7	1.12	0.91 - 1.33	0.07
171074	20.2	16.6 - 23.8	1.2	113	92 - 134	7	1.13	0.92 - 1.34	0.07
171075	20.2	16.6 - 23.8	1.2	113	92 - 134	7	1.13	0.92 - 1.34	0.07
172101	19.9	16.3 - 23.5	1.2	111	90 - 132	7	1.11	0.90 - 1.32	0.07
172104	19.9	16.3 - 23.5	1.2	111	90 - 132	7	1.11	0.90 - 1.32	0.07
173499	20.1	16.5 - 23.7	1.2	112	91 - 133	7	1.12	0.91 - 1.33	0.07
173501	20.1	16.5 - 23.7	1.2	112	91 - 133	7	1.12	0.91 - 1.33	0.07
175650	20.9	17.0 - 24.8	1.3	117	96 - 138	7	1.17	0.96 - 1.38	0.07
176136	21.0	17.1 - 24.9	1.3	117	96 - 138	7	1.17	0.96 - 1.38	0.07
176138	21.0	17.1 - 24.9	1.3	117	96 - 138	7	1.17	0.96 - 1.38	0.07
176469	20.4	16.8 - 24.0	1.2	114	93 - 135	7	1.14	0.93 - 1.35	0.07

LOT	Roche µmol/L			Roche µg/dL			Roche mg/L		
	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD
	Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE
	Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP
	Värde/Værdi	Område/Område	1DP/1SD	Värde/Værdi	Område/Område	1DP/1SD	Värde/Værdi	Område/Område	1DP/1SD
	Homota/Hodnota	Rozmezí/Rozsah	1SD/1SD	Homota/Hodnota	Rozmezí/Rozsah	1SD/1SD	Homota/Hodnota	Rozmezí/Rozsah	1SD/1SD
	Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD
Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	

### Precinorm U plus, No. 12149435 122

168539	19.4	15.8 - 23.0	1.2	108	90 - 126	6	1.08	0.90 - 1.26	0.06
172222	18.6	15.3 - 21.9	1.1	104	86 - 122	6	1.04	0.86 - 1.22	0.06
175487	18.5	15.2 - 21.8	1.1	103	85 - 121	6	1.03	0.85 - 1.21	0.06

### Precipath U, No. 10171760 122, 10171778 122, 10651265 922 (Q.C.S.)

167436	28.9	23.8 - 34.0	1.7	162	132 - 192	10	1.62	1.32 - 1.92	0.10
168009	30.8	25.4 - 36.2	1.8	172	142 - 202	10	1.72	1.42 - 2.02	0.10
168309	30.4	25.0 - 35.8	1.8	170	140 - 200	10	1.70	1.40 - 2.00	0.10
169839	30.0	24.6 - 35.4	1.8	168	138 - 198	10	1.68	1.38 - 1.98	0.10
169840	30.0	24.6 - 35.4	1.8	168	138 - 198	10	1.68	1.38 - 1.98	0.10
171368	31.1	25.4 - 36.8	1.9	174	144 - 204	10	1.74	1.44 - 2.04	0.10
173043	30.3	24.9 - 35.7	1.8	169	139 - 199	10	1.69	1.39 - 1.99	0.10
174531	30.5	25.1 - 35.9	1.8	170	140 - 200	10	1.70	1.40 - 2.00	0.10
175332	31.3	25.6 - 37.0	1.9	175	142 - 208	11	1.75	1.42 - 2.08	0.11
175333	31.3	25.6 - 37.0	1.9	175	142 - 208	11	1.75	1.42 - 2.08	0.11
176287	30.6	25.2 - 36.0	1.8	171	141 - 201	10	1.71	1.41 - 2.01	0.10

### Precipath U plus, No. 12149443 122

166448	31.2	25.5 - 36.9	1.9	174	144 - 204	10	1.74	1.44 - 2.04	0.10
171737	27.2	22.4 - 32.0	1.6	152	125 - 179	9	1.52	1.25 - 1.79	0.09
176975	29.6	24.2 - 35.0	1.8	165	135 - 195	10	1.65	1.35 - 1.95	0.10

LOT	µmol/L			µg/dL			mg/L		
	Sollwert	Bereich BÄK		Sollwert	Bereich BÄK		Sollwert	Bereich BÄK	

**Precinorm U, No. 10171735 122, 10171743 122, 10651257 922 (Q.C.S.)**

167326	20.7	17.7	-	23.7	116	101	-	131	1.16	1.01	-	1.31
167486	20.5	17.5	-	23.5	115	100	-	130	1.15	1.00	-	1.30
167969	21.9	18.9	-	24.9	122	104	-	140	1.22	1.04	-	1.40
167972	21.9	18.9	-	24.9	122	104	-	140	1.22	1.04	-	1.40
168314	20.6	17.6	-	23.6	115	100	-	130	1.15	1.00	-	1.30
168317	20.6	17.6	-	23.6	115	100	-	130	1.15	1.00	-	1.30
170571	20.0	17.3	-	22.7	112	97	-	127	1.12	0.97	-	1.27
170574	20.0	17.3	-	22.7	112	97	-	127	1.12	0.97	-	1.27
171074	20.2	17.5	-	22.9	113	98	-	128	1.13	0.98	-	1.28
171075	20.2	17.5	-	22.9	113	98	-	128	1.13	0.98	-	1.28
172101	19.9	17.2	-	22.6	111	96	-	126	1.11	0.96	-	1.26
172104	19.9	17.2	-	22.6	111	96	-	126	1.11	0.96	-	1.26
173499	20.1	17.4	-	22.8	112	97	-	127	1.12	0.97	-	1.27
173501	20.1	17.4	-	22.8	112	97	-	127	1.12	0.97	-	1.27
175650	20.9	17.9	-	23.9	117	102	-	132	1.17	1.02	-	1.32
176136	21.0	18.0	-	24.0	117	102	-	132	1.17	1.02	-	1.32
176138	21.0	18.0	-	24.0	117	102	-	132	1.17	1.02	-	1.32
176469	20.4	17.4	-	23.4	114	99	-	129	1.14	0.99	-	1.29

**Precinorm U plus, No. 12149435 122**

168539	19.4	16.7	-	22.1	108	93	-	123	1.08	0.93	-	1.23
172222	18.6	15.9	-	21.3	104	89	-	119	1.04	0.89	-	1.19
175487	18.5	15.8	-	21.2	103	88	-	118	1.03	0.88	-	1.18

**Precipath U, No. 10171760 122, 10171778 122, 10651265 922 (Q.C.S.)**

167436	28.9	25.0	-	32.8	162	138	-	186	1.62	1.38	-	1.86
168009	30.8	26.6	-	35.0	172	148	-	196	1.72	1.48	-	1.96
168309	30.4	26.2	-	34.6	170	146	-	194	1.70	1.46	-	1.94
169839	30.0	25.8	-	34.2	168	144	-	192	1.68	1.44	-	1.92
169840	30.0	25.8	-	34.2	168	144	-	192	1.68	1.44	-	1.92
171368	31.1	26.6	-	35.6	174	150	-	198	1.74	1.50	-	1.98
173043	30.3	26.1	-	34.5	169	145	-	193	1.69	1.45	-	1.93
174531	30.5	26.3	-	34.7	170	146	-	194	1.70	1.46	-	1.94
175332	31.3	26.8	-	35.8	175	151	-	199	1.75	1.51	-	1.99
175333	31.3	26.8	-	35.8	175	151	-	199	1.75	1.51	-	1.99
176287	30.6	26.4	-	34.8	171	147	-	195	1.71	1.47	-	1.95

**Precipath U plus, No. 12149443 122**

166448	31.2	26.7	-	35.7	174	150	-	198	1.74	1.50	-	1.98
171737	27.2	23.3	-	31.1	152	131	-	173	1.52	1.31	-	1.73
176975	29.6	25.4	-	33.8	165	141	-	189	1.65	1.41	-	1.89